

# THETA

The outstanding quality of the THETA re-usable product range is immediately apparent; high grade stainless steel, reassuring weight, excellent balance and superb finish. We strongly believe THETA is not only the best reusable product available worldwide but is also far superior to the majority of re-usable brands. For THETA products please refer to the product codes in blue.



# ALPHA

The disposable ALPHA range is as versatile as you need it to be, we will create any bespoke size and style upon request. The range shown in this catalogue is a guide to the most popular products we currently provide, many are regular stock items for fast delivery. Durable, high quality Stainless Steel shafts, reassuringly weighty and balanced Aluminium handles. For ALPHA products please refer to the product codes in gold.



## POWERED SUCTION

**Mercedes**  
Three holes spaced evenly around the tip.  
The most popular tip style.



ALPHA RANGE	Diameter (mm)	Length (cm)
ALPHA 210MM	2	10
ALPHA 315MM	3	15
ALPHA 325MM	3	25
ALPHA 430MM	4	30
ALPHA 530MM	5	30

THETA RANGE	Diameter (mm)	Length (cm)
THETA 210MM	2	10
THETA 315MM	3	15
THETA 325MM	3	25
THETA 430MM	4	30
THETA 530MM	5	30

**30 Hole harvest cannula**  
Harvest fat in smaller parcels, reduces blocking.



ALPHA RANGE	Diameter (mm)	Length (cm)
ALPHA 325/30hole	3	25

THETA RANGE	Diameter (mm)	Length (cm)
THETA 325/30hole	3	25

**Latern**  
Chinese lantern shaped tip, useful for fibrous areas and popular in Gynecomastia procedures.



ALPHA RANGE	Diameter (mm)	Length (cm)
ALPHA 210TSN	2	25
ALPHA 320TSN	3	30
ALPHA 430TSN	4	30

THETA RANGE	Diameter (mm)	Length (cm)
THETA 210TSN	2	25
THETA 320TSN	3	30
THETA 430TSN	4	30

**Submental**  
One hole, flat backed for areas closer to the skin



ALPHA RANGE	Diameter (mm)	Length (cm)
ALPHA 210SN	2	10
ALPHA 315SN	3	15

THETA RANGE	Diameter (mm)	Length (cm)
THETA 210SN	2	10
THETA 315SN	3	15

ALPHA RANGE	Diameter (mm)	Length (cm)
ALPHA 215OHN	2	15
ALPHA 320OHN	3	20
ALPHA 430OHN	4	30

THETA RANGE	Diameter (mm)	Length (cm)
THETA 215OHN	2	15
THETA 320OHN	3	20
THETA 430OHN	4	30



# FLUID INFUSION

## Infusion needles

High grade stainless steel. 6 spaced holes in a spiral around the tip. For fluid infusion using leur lock syringes, handles or giving set



ALPHA RANGE	Diameter (mm)	Length (cm)
ALPHA IMF/M 215	2	15
ALPHA IMF/M 223	2	23
ALPHA IMF/M 230	2	30
ALPHA IMF/M 315	3	15
ALPHA IMF/M 323	3	23
ALPHA IMF/M 330	3	30

THETA RANGE	Diameter (mm)	Length (cm)
THETA IMF/M 215	2	15
THETA IMF/M 223	2	23
THETA IMF/M 230	2	30
THETA IMF/M 315	3	15
THETA IMF/M 323	3	23
THETA IMF/M 330	3	30

## Curved



ALPHA RANGE	Diameter (mm)	Length (cm)
ALPHA IMF/M1.623CRV	1.6	23

THETA RANGE	Diameter (mm)	Length (cm)
THETA IMF/M1.623CRV	1.6	23

# NICRO HARVEST

**Open harvest**  
Guarded fully open tip, designed specifically for tissue harvest.



ALPHA RANGE	Diameter (mm)	Length (cm)
ALPHA NIC315/HARV	3	15
ALPHA NIC322/HARV	3	22

THETA RANGE	Diameter (mm)	Length (cm)
THETA NIC315/HARV	3	15
THETA NIC322/HARV	3	22

**30 Hole harvest cannula**  
Harvest fat in Micro droplets for fine graft re-injection.



ALPHA RANGE	Diameter (mm)	Length (cm)
ALPHA 2.515L/30hole	2.5	15
ALPHA 315L/30hole	3	15

THETA RANGE	Diameter (mm)	Length (cm)
THETA 2.515L/30hole	2.5	15
THETA 315L/30hole	3	15

**Submental**  
One hole, flat backed for areas closer to the skin.



ALPHA RANGE	Diameter (mm)	Length (cm)
ALPHA 322D/L	3	22

THETA RANGE	Diameter (mm)	Length (cm)
THETA 322D/L	3	22

**Mercedes tip**  
Three holes spaced evenly around the tip. The most popular tip style.



ALPHA RANGE	Diameter (mm)	Length (cm)
ALPHA 215DL	2	15
ALPHA 315DL	3	15
ALPHA 322DL	3	22

THETA RANGE	Diameter (mm)	Length (cm)
THETA 215DL	2	15
THETA 315DL	3	15
THETA 322DL	3	22

# TOOMEY SYRINGE SUCTION

## Mercedes tip

Three holes spaced evenly around the tip. The most popular tip



ALPHA RANGE	Diameter (mm)	Length (cm)
ALPHA 210NTS	2	10
ALPHA 315NTS	3	15
ALPHA 318NTS	3	18
ALPHA 325NTS	3	25
ALPHA 430NTS	4	30

THETA RANGE	Diameter (mm)	Length (cm)
THETA 210NTS	2	10
THETA 315NTS	3	15
THETA 318NTS	3	18
THETA 325NTS	3	25
THETA 430NTS	4	30

## 30 Hole harvest cannula

Harvest fat in smaller parcels, reduces blocking.



ALPHA RANGE	Diameter (mm)	Length (cm)
ALPHA2.515/30holeS	2.5	15
ALPHA 320/30holeTS	3	20

THETA RANGE	Diameter (mm)	Length (cm)
THETA2.515/30holeTS	2.5	15
THETA320/30holeTS	3	20

## Submental

One hole, flat backed for areas closer to the skin.



ALPHA RANGE	Diameter (mm)	Length (cm)
ALPHA 210SNTS	2	10
ALPHA 315SNTS	3	15

THETA RANGE	Diameter (mm)	Length (cm)
THETA 210SNTS	2	10
THETA 315SNTS	3	15

## Open harvest

Guarded fully open tip, designed specifically for tissue harvest.



ALPHA RANGE	Diameter (mm)	Length (cm)
ALPHA 318HARV/TS	3	18

THETA RANGE	Diameter (mm)	Length (cm)
THETA 318HARV/TS	3	18

## NICRO INFILTRATION CANNULAE

High grade Stainless Steel for re-injection using Luer lock syringes



ALPHA RANGE	Diameter (mm)	Length (cm)
ALPHA NIC16/A7	16G	7
ALPHA NIC16/A9	16G	9
ALPHA NIC16/B7	16G	7
ALPHA NIC16/B9	16G	9
ALPHA NIC16/B12	16G	12
ALPHA NIC16/C7	16G	7
ALPHA NIC16/C9	16G	9
ALPHA NIC18/A9	18G	9
ALPHA NIC18/B9	18G	9

THETA RANGE	Diameter (mm)	Length (cm)
THETA NIC16/A7	16G	7
THETA NIC16/A9	16G	9
THETA NIC16/B7	16G	7
THETA NIC16/B9	16G	9
THETA NIC16/B12	16G	12
THETA NIC16/C7	16G	7
THETA NIC16/C9	16G	9
THETA NIC18/A9	18G	9
THETA NIC18/B9	18G	9

### Curved Cannula



ALPHA RANGE	Diameter (mm)	Length (cm)
ALPHA CURV/18B7	18G	7
ALPHA CURV/16B12	16G	12
ALPHA CURV/16B15	16G	15

THETA RANGE	Diameter (mm)	Length (cm)
THETA CURV/18B7	18G	7
THETA CURV/16B12	16G	12
THETA CURV/16B15	16G	15